

Required Components



BOM

- 1. Grip Handle
- 2. Grip insert
- 3. Insert
- 4. Spoon attachment
- 5. Bearing
- 6. Fishing weights
- 7. Utensil
- 8. Computer case screw

Required Tools

- Tin Snips
- Superglue
- Hot Glue

Required Personal Protective Equipment (PPE)

Eye protection



Assembly Instructions

Step 1

Use superglue to glue the bearing into the grip handle. Be careful to not get any into the bearing.



Step 2
Use superglue to attach the insert into the bearing. Be careful not to get any into the bearing.





Step 3
Cut the utensil at the desired length and hot glue it into the slot on the spoon attachment.



Step 4
Using the set screw, attach the spoon section to the handle section.





Step 5

If the user wants a weighted utensil, add the desired amount of weights to the hollow handle section. A tissue can optionally be used to prevent the weights from sliding around in the handle.



Step 6
Use superglue to attach the grip insert to the end of the handle. Optionally, a soldering iron can be used to melt the two pieces together.

