

Utensil Holder

ASSEMBLY MANUAL

MATERIALS



Utensil holder



3/4" Velcro




Scissors



Hot glue




Any utensil



01. Measure and cut a piece of Velcro to fit in the shorter strap indent.



02. Apply hot glue.



03. Stick the Velcro with the **hook (rough)** side exposed.



04. Apply hot glue in the longer strap indent.




05. Stick the Velcro with the **loop (smooth)** side exposed.



06. Hold the utensil holder and secure the Velcro to the user's comfort.

Utensil Holder

ASSEMBLY MANUAL



07. Cut the Velcro strap to length.



08. Insert a utensil.

