

# DRAG

## Maker Guide

© 2024 PrintLab International Ltd. Licensed under [CC BY SA 4.0](#)




v2. 22nd February 2024 – screw holes in left and centre plates  
updated to clearance holes



DRAG is an assistive device that enables users to write or draw without having to form a tight closed grip with their fingers. Simply rest your hand on the ergonomic mouse-shaped body, then press and DRAG. The device is suited for people with arthritis or other hand mobility issues that cause pain or discomfort when writing or drawing. Its symmetrical design also caters to the needs of both left and right-handed users.

Video of device in action: <https://bit.ly/47RFs2M>

Video of design process and rationale: <https://bit.ly/47RFOGE>



≈Weight (excluding pen): 78g

# MATERIALS


Description	Qty Required	≈ Cost (USD)	Example Options
M3 X 20mm Machine Screws	3	\$1.68 (pack of 6)	<a href="https://bit.ly/4bbDLzT">https://bit.ly/4bbDLzT</a>
6mm Diameter Mouse Skates <small>*Optional but recommended for smoother operation</small>	6	\$9.99 (pack of 20)	<a href="https://bit.ly/3HEc5Gy">https://bit.ly/3HEc5Gy</a>
3D Printed Left Plate	1	\$0.86 (PLA material cost)	-
3D Printed Centre Plate	1	\$0.27 (PLA material cost)	-
3D Printed Right Plate	1	\$0.85 (PLA material cost)	-
Pen (up to 12mm diameter)	1	-	-

---

## Tools Required


3D Printer with PLA filament

Screwdriver (Type dependent on screws used)



# 3D PRINTING

All 3 files may be 3D printed on the same print bed.  
Please follow the below guidance on orientation and  
print settings.



## Settings


Technology	FFF
Material	PLA
Nozzle Diameter	0.4mm
Layer Height	0.2mm
Infill	15%
Support	None

## Estimated Stats

Print Time	4 hours
Mass	76g
Cost	\$2 (USD)

# ASSEMBLY

---



## Parts List



- A. 1 x 3D printed right plate
- B. 1 x 3D printed centre plate
- C. 1 x 3D printed left plate
- D. 3 x M3x20 machine screws
- E. 6 x 6mm diameter adhesive mouse skates (optional)
- F. 1 x pen

## Tools Required

Screwdriver  
(Type dependent on screws used)


# 01.

Press the 3 plates together, ensuring the shapes and holes align with each other.  
Then connect them with 2 screws.




# 02.

(Optional Step). Place a series of adhesive mouse skates onto the underside of the left plate and right plate. Position them on the inner edges where the device will meet the paper or surface it is being used on.



# 03.

Place a pen into the groove and tighten a screw through the left plate and into the right plate. When you begin to feel the groove tighten around the pen, slide then pen so the tip is just above the writing/drawing surface when placed on its base. Tighten the screw further until the pen is secure in position.




# USAGE



DRAG can be used for both writing and drawing. Simply rest your hand on the mouse-shaped body and place 2 fingers either side of the pen. To activate the pen, gently press your 2 fingers down and drag to write or draw. You may also experiment with different pen positions to suit your needs. To do this, loosen the screw and slide the pen to a different position, before retightening.



# EVOLUTION



The DRAG design files have been released under a Creative Commons [CC BY SA 4.0](https://creativecommons.org/licenses/by-sa/4.0/) License, whereby we encourage people to modify, adapt and improve the design as they see fit. We are particularly interested to see how the design might be refined to be more customisable for different users, as well as how the parts could be consolidated. If you have any feedback or suggestions, please email them to [hello@weareprintlab.com](mailto:hello@weareprintlab.com).

[www.weareprintlab.com](http://www.weareprintlab.com)