


Spoon Stabilizer Support

USER QUICK GUIDE

Introduction

A device for those who have difficulties with tremors when eating. It uses both passive stabilization and weights to dampen the effect of hand tremors.

Features



Usage

Attach either the spoon or fork attachment to the handle using the set screw. The device can then be used as a normal spoon or fork.

Compatibility

The device is compatible with any spoon or fork that has a cross section of 9x2.5mm or less.

Specifications

Handle has a length of 90mm, and a diameter of 30mm. The entire device has a weight of 60 g.

Cleaning

The device can be wiped clean with a damp cloth. The section with the fork/spoon can be detached for easier cleaning.



© 2023 by [Neil Squire](#).

This work is licensed under the CC BY SA 4.0 License: <http://creativecommons.org/licenses/by-sa/4.0>

Files available at <https://makersmakingchange.com/project/spoon-stabilizer-support/>